

VARIABEL WINKELSTABILE DISTALE RADIUSPLATTE

Variable angle-stable distal radial plate

- Titan

MINIFRAGMENT 2,7 mm

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Nearly all fractures of the distal radius can be managed surgically with the 2.7 mm volar fixed-angle radius plate for internal fixation. The anatomical shape of the plate in combination with the variable angular stability of the two distal rows of holes enables stable management of all extra-articular and nearly all intra-articular fractures (C1-3 in the AO classification) without additional measures. Particularly in the case of intra-articular comminuted fractures, the variable angular stability provides excellent conditions for adjusting the screw position to the individual fracturesituation. The small holes incorporated in the plate can be used for preliminary introduction of K-wires, which makes the intraoperative procedure considerably easier. Preliminary K-wire fixation allows the position of the plate relative to the radius and the reduction result to be assessed easily under image intensifier control and corrected if necessary. All of the required plate lengths are available, each in two different widths at the distal end (23 mm and 27 mm).

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Treatment with 2.7 mm variable fixed-angle 7/2-hole radius plate.

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Positioning splint, where applicable, easily removable wrist orthosis with Velcro fastening. Start independent exercises without range of motion restrictions from the first postoperative day. Physiotherapy is required in only a few cases.

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The operation is performed in supine position. The arm is placed on an arm table. Closed reduction is carried out first, if necessary under image intensifier control. In A3 fractures and the majority of C2 fractures, anatomical reduction is largely achieved when ligamentotaxis is maintained.

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The operation is performed without a tourniquet and meticulous haemostasis is obtained with bipolar coagulation. Short (approx. 3-5 cm) incision over the flexor surface of the wrist. Exposure of the flexor carpi radialis tendon, which is retracted slightly toward the ulnar side.

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
The forearm fascia is divided immediately below or just to the radial side of the tendon. The flexor carpi radialis tendon and the finger flexor tendons / muscle bellies beneath it are held aside (using a self-retaining retractor if necessary). The median nerve is not exposed and the carpal tunnel is not touched (exception: pre-existing CTS). Division of pronator quadratus with diathermy on the radial side, and the muscle belly is pushed away from the flexor surface of the radius.

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The fracture line on the flexor surface is exposed cleanly by sparing removal of periosteum. This is followed by definitive anatomical reduction.

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
A suitable plate is selected (the two-hole plate is usually sufficient, 23 mm or 27 mm wide, depending on the width of the radius).

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The plate is fixed to the articular fragment and to the shaft with one or more Kirschner wires through the holes provided in the plate. This achieves preliminary fixation of the plate to the radius and of the reduction (image intensifier control). The plate position is corrected if necessary with further preliminary Kirschner wire fixation and subsequent image intensifier control.

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During drilling, the opposite cortex is deliberately drilled cautiously with a finger on the dorsal side of the radius. The length is measured through the fixed-angle drill bushing or after using the variable fixed-angle drill bushing with corresponding distance sleeve on the length gauge, palpating the tip of the gauge on the dorsal side of the radius beneath the skin so that the length of the screws can be established exactly.

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
The screws introduced distally should just reach the dorsal cortex but not push away smaller cortical fragments (which are common when there is a dorsal zone of comminution). It is therefore beneficial, as mentioned above, to drill the thin dorsal cortex carefully first. In this technique, at least 3 or 4 screws are introduced in succession into the articular block and the radial screw in particular should diverge into the radial styloid process depending on the inclination. In the case of intra-articular fractures, a screw position in the fracture gap should be avoided.

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
Important: when tightening the distal fixed-angle screws, counterpressure should be exerted with a finger from the dorsal side in order to press the bone against the plate. That way, the anatomical shape of the plate is used for exact reduction.

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Final image intensifier control and documentation in two planes.

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After X-ray documentation and irrigation, insertion of a vacuum drain and skin closure.

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Treated with 5/2-hole variable
fixed-angle radius plate

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Extension fracture with
impression of the articular
surface, type C3

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Treated with 7/2-hole variable
fixed-angle radius plate

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
The postoperative CT shows an anatomically reconstructed articular surface without a step, with the screws in correct, strictly subchondral (extra-articular!) position.

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Early functional follow-up treatment

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Range of motion on the 3rd postoperative day, early functional follow-up treatment.
No immobilisation!

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Extension fracture type C2 with depressed ulnar marginal fragment

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Reduction of the ulnar marginal fragment with Kirschner wire from the dorsal side ("joy stick"), at the same time preliminary fixation of this fragment with a Kirschner wire from the volar aspect through the plate.

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"Avoiding" the intra-articular fracture gap by exploiting the variable angular stability. Treated with variable fixed-angle radius plate 7/3-hole.

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Juxta-articular forearm fracture (70-year old patient, radius flexion fracture).
Treated with 5/3-hole plate on the radius and 2.0 mm distal fixed-angle ulnar plate.

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Comminuted fracture (C3) 38-year old patient.
Healed appearance 6 months postoperatively, metal removal planned.

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Juxta-articular forearm fracture with C2 fracture of the radius, pre-existing bayonet deformity after old radius fracture healed in deformity (78-year old patient). Corrective fixed-angle internal fixation of the radius with 7/3-hole plate and fixed-angle internal fixation of the fracture of the ulnar head with 2.0 mm distal ulnar plate.

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

Distal forearm comminuted fracture, high-grade osteoporosis (85-year old patient).
Long bridging, volar supporting, fixed-angle internal fixation of the radius with 5/8-hole plate, postoperative immobilisation in a cast.

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Order additionally

Winkelstabile distale Radiusplatte schmal

Fixed-angle stable distal radial plate small

Bestell-Nr. Titan

Code N° titanium

5.9823.52	rechts / right
5.9833.52	links / left
5.9823.53	rechts / right
5.9833.53	links / left
5.9823.54	rechts / right
5.9833.54	links / left
5.9823.56	rechts / right
5.9833.56	links / left
5.9823.58	rechts / right
5.9833.58	links / left
5.9823.510	rechts / right
5.9833.510	links / left
5.9823.513	rechts / right
5.9833.513	links / left

Länge

Length

43 mm
43 mm
51 mm
51 mm
60 mm
60 mm
77 mm
77 mm
94 mm
94 mm
111 mm
111 mm
136 mm
136 mm

Lochanzahl

N° of hole

2
2
3
3
4
4
6
6
8
8
10
10
13
13

Winkelstabile distale Radiusplatte breit

Fixed-angle stable distal radial plate

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Cortical screws, diameter 2,7 mm (self-tapping)

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
Cortical screws, diameter 2,7 mm, conical head thread (self-tapping)

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19.115.00

1

Locating sieve (for instruments and implants)

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2.940.25

1

Hexagonal screw-driver, with handle, for screws diameter 2,7 mm

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2.904.03

1

Spiral drill (for quick coupler) diameter 2,0 mm, length 105 mm

Schraubenhaltepinzette

2.954.01

1

Screw-holding forceps

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2.953.65

1


Screw gauge

Kirschnerdraht Trokar, Ende rund, D 1,6 mm

6.031.16


5

Kirschner-wire with trokar point, rounded end, D 1,6 mm

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No part of this publication may be reproduced, stored in a retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying, recording or otherwise without the prior written permission of the publisher.

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© Copyright Königsee Implantate GmbH
OSTEOSYNTHESE GmbH
Ltd. OA. Dr. med. Klaus Jürgen Maier
ZUR OSTEOSYNTHESE GmbH**